
Industrial Division

VR2 Training

1 | September 30, 2016 | Proprietary © 2014 Curtiss-Wright

Curtiss-Wright

Remote Controller nVR2

VR2 system

- Modular wheelchair control system
- VR2 or R-net type motor and battery connections
- Drive plus 1 or 2 actuators
- Easy installation on all types of wheelchairs
- Service-friendly

VR2 is a multi-module control system designed specifically to cater for the needs of entry-level to mid-range wheelchairs. The VR2 offers a high degree of flexibility, including options for the control of up to two seating actuators. In addition, versions of the Power Module with R-net connectors are available, so allowing simple upgrade paths for wheelchair models or users which may need to transition from a VR2 based system.

The Power Module is available with R-net type motor and battery connectors offering simple upgrade path on a wheelchair that requires more complex seating controls or specialty inputs available with R-net.

Joystick Modules

- Ergonomic and stylish design
- Drive only, seating and lighting options
- Simple and clear user controls
- Service friendly
- Thin and flexible cable

There are multiple VR2 Joystick Modules available - from the drive-only version to options that can control seating and lighting functions - each being extremely ergonomic and stylish. The user controls are intuitive and easy to understand, making it simple for new wheelchair users to quickly become confident. All versions have been designed with servicing in mind; for example, the keypads and cable are field replaceable.

Dual Attendant Module

- Push button exchange of control between user and attendant
- Dual Attendant control of max. speed & actuators
- Works with all Power Modules as standard

The VR2 Dual Attendant Control system consists of a Joystick Module, Power Module and Dual Attendant Module. This allows the drive and seat actuator functions of the wheelchair to be controlled by either the occupant or by an attendant from another location on the wheelchair. Control can easily be exchanged via a push button and the attendant can limit the speed of the wheelchair. The attendant also has the option to adjust seat actuators when connected to a suitable VR2 power module.

VR2 Product Family

Power Modules

- Choice of power ratings up to 90A max
- Drive plus 1 or 2 actuators
- VR2 or R-net type motor and battery connections as standard
- Mounting dimensions as per R-net
- Versatile Inhibits

To offer the optimum combination of power and price the VR2 Power Module is available in 60A, 70A and 90A models. To add further flexibility to the system each Power Module can be populated to drive one or two seat actuators and is configured to recognize a Lighting and Dual Attendant Modules as standard. Each Power Module is available with R-net type battery and motor connectors and has the same mounting dimension as an R-net Power Module. This offers the flexibility to standardize motor and battery harnessing where R-net is a configuration option for the particular wheelchair. Advances in software design have led to a gradient performance never previously available in cost effective systems. If the wheelchair is driven up a hill there is practically zero roll-back between the time the joystick is released and the brakes are applied.

A high-power connector, primarily intended for an on-board charger, is fitted to the Power Module as standard. The 12Amps rating of this connector means it can be safely used for charging large capacity batteries. A universal inhibit pin is included, which can be set up with any type of charger inhibit output. A second inhibit connector is also provided and can also be used wherever there is a need to limit or prevent drive or seat movements.

Lighting Module

- Separate module to control lighting functions
- Works with all Power Module as standard
- StVZO lighting control
- Small foot print
- Modular system that is easy to install


The VR2 lighting system consists of a Joystick Module, Power Module and Lighting Module. This configuration offers a cost effective elegant way to control lights, turn indicators and hazard lights. The system provides constant brightness lamp control and turn indicator fault detection in compliance with the German StVZO standard, which is mandatory for wheelchairs in that market place.

The small size and modular nature of the lighting system allows for unobtrusive installations anywhere on the wheelchair. The simple push fit connection system makes wiring extremely quick and easy. With self protecting outputs there is no need to fit external fuses or breakers, so wiring harnesses are kept simple.

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




The table below gives the minimum recommended wire sizes defined in ISO7176: 2008.

Controller Current Limit (A)	Battery wire size (mm ²)		Motor Wiring size (mm ²)	
	For Length <1000mm	For length 1000mm - 1500mm	For Length <1000mm	For length 1000mm - 1500mm
40	3.0	4.0	2.5	3.0
50	4.0	4.0	3.0	4.0
60	6.0	6.0	3.0	4.0
80	6.0	8.0	4.0	5.0
100	6.0	No guidance	6.0	6.0

Connector Kits	PGDT Reference	Inconnect Reference
Motor	D50750	IPG-7401-PS
Battery	D50751	IPG-7401-PS
Motor (90)	D49713 PG80-M	IPG-5401
Battery (90)	D49712 PG80-B	IPG-5402
On Board Charger	D50752	IPG-7301-PS
Inhibit 2	D50753	IPG-8201-PS
Actuator	D50754	IPG-8202-PS

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With the control system switched off, check that the joystick is not bent or damaged and that it returns to the center when you push and release it. If there is a problem do not continue with the safety checks and contact Your service agent.

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每週檢查 (Weekly Check)

電磁煞車 (Solenoid brake)

本測試需在水平地面進行，請確認輪椅至少須與附近障礙物保持1M以上的距離。
This test should be carried out on a level with at least one meter clear space around the wheelchair.

- 打開控制器。(Switch on the control system)
- 間隔1秒後，檢查電量表是否恆亮或緩慢閃爍。
(Check That the battery gauge remains on, or flashes slowly, after one second)
- 慢慢將搖桿向前推動，直到你聽到電磁煞車釋放的聲音。此時輪椅應可開始移動。
(Push the joystick slowly forwards until you hear the solenoid brakes operate. The chair may start to move)
- 立刻放開搖桿，你應該可以在短時間內聽到兩邊電磁煞車操作的聲音。
(Immediately release the joystick. You must be able to hear each solenoid brake operate within a few seconds)
- 依相同方法，分別測試搖桿向後、向左及向右推動的狀況。
(Repeat the test a further three times, pushing the joystick slowly backwards, left and right)

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Make sure that all connectors are securely mated.

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Check the condition of all cables and connectors for damage.

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Check the thin rubber gaiter or boot, around the base of the joystick shaft, for damage or splitting. Check visually, do not handle the gaiter.

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Make sure that all the components of the control system are securely mounted. Do not overtighten any securing screws.

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(Maximum Speed / Profile Indicator is Steady)

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The display will vary slightly depending on whether the control system is programmed to operate with drive profiles.

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(Maximum Speed / Profile Indicator Flashes)

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This indicates the speed of the wheelchair is being limited for safety reasons.

The exact reason will depend on the type of wheelchair, however, the most Common cause is that the seat is in the elevated position.

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(Maximum Speed / Profile Indicator Ripples Up and Down)

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This indicates the control system is locked.

輪椅控制系統解鎖步驟 (Unlocking the Wheelchair)

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Use the on/off button to switch the control system on. The maximum speed / profile indicator will be rippling up and down.

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Deflect the joystick forwards until the control system beeps.

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Deflect the joystick in reverse until the control system beeps.

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Release the joystick, there will be a long beep.

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The wheelchair is now unlocked.

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電量表穩定顯示不閃爍 (Battery Gauge is Steady)

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This indicates that all is well.

電量表緩慢閃爍 (Battery Gauge Flashes Slowly)

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The control system is functioning correctly, but you should charge the battery as soon as possible.

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The wheelchair batteries are being charged. You will not be able to drive the wheelchair until the charger is disconnected and you have switched the control system off and on again.

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(Battery Gauge Blinks Once Every 2.5 Seconds)

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The control system has “gone to sleep” because the wheelchair has not been driven for a period of time. The time period depends on the programming of the system. To re-start, switch the system off and on again.

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(Battery Gauge Flashes Rapidly even with the joystick released)

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This indicates a system trip, i.e. the controller has detected a problem somewhere in the wheelchair's electrical system.

Please follow the following procedure.

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Make sure that all connectors on the wheelchair and the control system are mated securely.

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If you can't find the problem, try using the self-help guide.

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
Switch on the control system again and try to drive the wheelchair. If the safety circuit operate again, switch off and do not try to use the wheelchair. Contact your service agent.

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1 Bar 	低電壓 Battery Need Charging	7 Bar + S 	操縱桿通訊不良 Communication Fault	Blink 	休眠狀態 Sleep Mode
2 Bar 	左馬達斷線 L.Motor Bad Connection	8 Bar + A 	電動缸故障 An Actuator Trip	Charge Step 	充電中 Charging
3 Bar 	左馬達短路 L.Motor Short Circuit			Ripple 	操縱桿未置中 Joystick Not Centralized
4 Bar 	右馬達斷線 R.Motor Bad Connection				
5 Bar 	右馬達短路 R.Motor Short Circuit				
6 Bar 	充電限速中 Speed Inhibit, Could be Charging Inhibit				
7 Bar 	操縱桿可能故障 Joystick Fault				
8 Bar 	控制器可能故障 Controller Fault				
9 Bar 	電磁煞車故障 S.Brake Bad Connection				
10 Bar 	電池電壓過高 Excessive Voltage, Battery Poor Connection				

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Adjusts the value for forward and reverse acceleration of the Wheelchair.

Adjustable in steps of 1 from 0 to 100.

A higher value gives faster acceleration.

This programmed value of acceleration occurs when speed setting 5 is selected.

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Adjusts the value for forward and reverse deceleration (or braking) of the wheelchair.

Adjustable in steps of 1 from 0 to 100.

A higher value gives faster deceleration.

This programmed value of deceleration occurs when speed setting 5 is selected.

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Adjusts the value for turning acceleration of the wheelchair,
from 0 to 100 in steps of 1.

A higher value gives faster acceleration. This programmed value of
acceleration occurs when speed setting 5 is selected.

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Adjusts the value for turning deceleration (or braking) of the wheelchair.
Adjustable in steps of 1 from 0 to 100.

A higher value gives faster deceleration. This programmed value of
deceleration occurs when the Joystick Module has speed setting 5 selected.

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Adjusts the minimum and maximum values for forward speed of the
wheelchair.

Adjustable in steps of 1% from 0 to 100%.

A higher value gives a faster speed. The minimum value occurs when speed
setting 1 is selected, and the maximum value occurs at speed setting 5.

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Adjusts the minimum and maximum values for reverse speed of the wheelchair.

Adjustable in steps of 1% from 0 to 100%.

A higher value gives a faster speed. The minimum value occurs when the speed setting 1 is selected, and the maximum value occurs at speed setting 5.

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Adjusts the minimum and maximum values for the turning speed of the wheelchair. Adjustable in steps of 1% from 0 to 100%.

A higher value gives a faster speed.

If Drive Profile 0 is selected then the minimum value occurs when speed setting 1 is selected, and the maximum value occurs when speed setting 5 is selected.

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Sets the power of the wheelchair.

Adjustable in steps of 1% from 0 to 100%.

Power is the ability of a wheelchair to climb a hill or overcome an obstacle.

If it is set to 100% then the wheelchair will provide full power.

Values below 100% will result in reduced power.

A typical use is to minimize damage to doorways or furniture if the wheelchair is being used indoors.

The values can be set independently between drive profiles, meaning separate indoor and outdoor profiles can be defined.

Speed

Number of Drive Profiles (Profiles ?) 行駛情境

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The controller can operate with single or multiple drive profiles.

A drive profile is a collection of programmable parameters comprising of Acceleration, Deceleration, Turn Acceleration, Turn Deceleration, Forward Speed, Reverse Speed and Turning Speed.

The number of drive profiles is determined by the programmable parameter, Number of Drive Profiles.

If Number of Drive Profiles is set to 0, then there is one setting for each of the parameters listed previously, and the control system's maximum speed setting can be changed with the maximum speed / profile increase and decrease buttons.

If the values of Number of Drive Profiles is 2 to 5, there is a corresponding number of drive profiles and each listed parameter can be individually set within a profile. The maximum speed /profile increase and decrease buttons are then used to switch between the available profiles.

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Actuator

Number of Actuators (Engineer Menu -> Actuators) 電動缸數量 設定使用電動缸的數量

Sets the number of actuators to be connected to the VR2 control system.

This value must be set correctly in order to get the most informative and

user-friendly method of actuator control.

It is adjustable between 0 - 2 actuators in steps of 1.

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Because of the wide variation in actuator motors and applications, the end-stop current threshold is programmable.

Actuator 1 Current Limit Timeout, Actuator 2 Current Limit Timeout 電動缸輸出達到停止力量後的延遲時間

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Actuator

Actuator Selection with Joystick (Act Selection? -> with J' stick)

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This parameter determines whether selection between different actuator

channels can be made with left/right movements of the joystick.

You can set between on and off.

On: Means that selection between actuator channels can be made with the joystick and the actuator buttons.

Off: Means that only the actuator buttons can be used to select between actuator channels.

The parameter is only applicable if the Number of Actuators parameter is set to greater than 1.

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This allows the use of an Audible bleep to tell you when the Actuator is stalled at its end stop.

This parameter can be set to On or Off.

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Operation

Sleep Time (minutes) (Sleep Timer ?) 休眠計時

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Sets the period of time before the control system will turn itself off if the wheelchair is not driven.

The time can be set between 0 and 30 minutes in steps of 1 minute.

If the time is set to 0 the system will never turn itself off.

Operation

Joystick Throw -> Forward (%) (Joystick Throw ? -> forward thrw)

Joystick Throw -> Reverse (%) (Joystick Throw ? -> reverse thrw)

Joystick Throw -> Left (%) (Joystick Throw ? -> left throw)

Joystick Throw -> Right (%) (Joystick Throw ? -> right throw)

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This allows you to program the control system so that full speed can be reached with a reduced joystick movement (throw).

This is particularly useful for wheelchair users with limited hand or arm movement.

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On Deflecting the joystick Forward will result in Reverse drive.

Off Deflecting the joystick Forward will result in Forward drive.

Left and Right deflection of the joystick remain unchanged.

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This parameter compensates for any mismatching of motors to ensure that the wheelchair drives directly forward when the control system's joystick is being pushed directly forward.

It is normally set to zero but may be varied from -9 to +9 in steps of 1.

If the chair is veering to the left, you should increase the setting.

If the chair veers to the right, decrease the setting.

If Swap Motors is set, this logic will be reversed.

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Switch Bleep Volume sets the volume of the audible feedback given whenever a button on the controller is operated.

Adjustable between 0 and 10 in steps of 1.

If the parameter is set to 0 then the Switch Bleep function is effectively switched OFF.

The higher the value, the louder the audible feedback.

Operation

Speed Adjustment while Driving (Speed Adjust ? -> in drive)

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This parameter sets whether the controller's speed/profile buttons are active while the wheelchair is being driven.

The parameter can be set to on or off.

On Means the buttons are active while the wheelchair is being driven, so the user can make maximum speed setting adjustments (or select a different drive profile) while actually moving.

Off Means the buttons are not active while the wheelchair is being driven, so the joystick must be released and the wheelchair at rest before maximum speed setting adjustments (or different drive profile selections) can be made.

Battery

Low Battery Flash Level (Engineer Menu ? -> Battery Flash ?)

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Low Battery Alarm (Engineer Menu ? -> Battery Alarm ?)

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

- Inhibit 2 Upper Level Threshold: 54%
- Inhibit 2 Middle Level Threshold: 34%
- Inhibit 2 Lower Level Threshold: 12%
- Inhibit 2 Speed Limit in Band 0: 100% 座位高低改變時改變
- Inhibit 2 Speed Limit in Band 1: 75% 行駛速度
- Inhibit 2 Speed Limit in Band 2: 50%
- Inhibit 2 Speed Limit in Band 3: 0%

VR2 ACTUATOR INHIBIT POSSIBLE SWITCH SET-UP

VR2 THRESHOLD - BAND RELATIONSHIP

Band 3	100%
Band 2	54%
Band 1	34%
Band 0	12%
	0%

Inhibit

Channel 1 Up

Inhibit in Inhibit 2 Band 0	Yes/No
Inhibit in Inhibit 2 Band 3	Yes/No
Inhibit in Inhibit 3 Band 0	Yes/No
Inhibit in Inhibit 3 Band 3	Yes/No

Channel 1 Down

Inhibit in Inhibit 2 Band 0	Yes/No
Inhibit in Inhibit 2 Band 3	Yes/No
Inhibit in Inhibit 3 Band 0	Yes/No
Inhibit in Inhibit 3 Band 3	Yes/No

Channel 2 Up

Inhibit in Inhibit 2 Band 0	No
Inhibit in Inhibit 2 Band 3	Yes
Inhibit in Inhibit 3 Band 0	No
Inhibit in Inhibit 3 Band 3	No

Channel 2 Down

Inhibit in Inhibit 2 Band 0	No
Inhibit in Inhibit 2 Band 3	Yes
Inhibit in Inhibit 3 Band 0	No
Inhibit in Inhibit 3 Band 3	No

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Disable Inhibit x in Drive 行駛間忽略Inhibit的變化

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General (PM90 only)

Front Wheel Drive Type , Front Wheel Drive Rate (Engineer Menu ? -> Front Drv Rate ?) , FWD Acceleration Limit , Advanced Electronic Stability (AES) AES Rate 前驅修正參數

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
Motor

Current Foldback Temperature (Deg C)

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
設定範圍: 1 ~ Current Limit Max

Motor

Current Foldback Threshold 電流保護門檻

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
Current Foldback Level

(Engineer Menu ? -> Timed Foldback ? -> level)

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Motor

Motor Compensation (Engineer Menu ? -> Compensation ?)

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Demand clipping is necessary to ensure good speed and direction controller of the wheelchair in all circumstances, including heavy loads and low battery conditions.

Off If set to Off it uses the standard controller algorithm

On If set to On it uses the R-net version of the demand clipping algorithm.

If some wheelchair application setting the parameter On can allow higher speeds to be maintained when making regular course adjustments.

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The Torque parameter boosts the current to the motors at low speed settings.

If the motor is stalled, for example, the wheelchair is stuck against an obstacle, such as a door threshold; then this will be automatically detected and the current to the motors will be increased, allowing the obstacle to be overcome.

Torque can be set between 0% and 100%

A value of 0% means the Torque parameter has no effect.

Higher values mean that more current will be permitted in the described stall conditions.

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This parameter allows the effects of hand tremor to be reduced. If the user has a condition that results in hand tremor, then increasing the value of Tremor Damping will reduce the effect of the tremor, making the wheelchair more controllable.

Tremor Damping can be set between 0% and 100%

A value of 0% means Tremor Damping has no effect.

Note, even at this value, there is inherent damping in the control system. Higher values apply a higher level of The higher Tremor Damping is set the slower joystick response will become.

Joystick Stationary Time (Mins) 搖桿靜置計時器

Joystick Stationary Range 搖桿靜置的容許誤差範圍

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Ripple

操縱桿未置中
Joystick Not Centralized

Trip Code	Trip Type	Description & Reference
0A00	-	Controller in Sleep Mode - 3.18
1320	-	Timed Foldback Active - 3.16
1D05	7	JS Time Exceeded - 3.7
1500	9	Solenoid Brake Trip - 3.9
1505	9	Left Solenoid Brake Trip - 3.9
1506	9	Right Solenoid Brake Trip - 3.9
1600	10	High Battery Voltage - 3.10
1E03	Charging	Refer to section - 3.6
1E04	6	Refer to section - 3.13
1E05	Charging	Refer to section - 3.14
2C00	1	Low Battery Voltage - 3.1
2C02	-	Low Battery Lock-out - 3.1
2F00	User	Refer to sections 3.7 & 3.11
3B00	2	Left Motor Disconnected - 3.2
3C00	4	Right Motor Disconnected - 3.4
3D00	3	Left Motor Wiring Trip - 3.3
3D01	3	Left Motor Wiring Trip - 3.3
3E00	3	Right Motor Wiring Trip - 3.3
3E01	3	Right Motor Wiring Trip - 3.3
4401	8	Control System Trip - 3.8

Trip Code	Trip Type	Description & Reference
5400	7+5 ¹	Communications Trip - 3.12
7A03	A Only ¹	Actuator Motor Wiring Trip - 3.15
7100	7	Joystick Trip - 3.7
7101	7	Joystick Trip - 3.7
7102	7	Joystick Trip - 3.7
7103	7	Joystick Trip - 3.7
7104	7	Joystick Trip - 3.7
7105	7	Joystick Trip - 3.7
7107	7+0 ¹	Refer to sections 3.19
7140	7+0 ¹	Attendant Module Trip 3.19
7141	7+0 ¹	Attendant Module Trip 3.19
7142	7+0 ¹	Attendant Module Trip 3.19
7143	7+0 ¹	Attendant Module Trip 3.19
7145	7+0 ¹	Attendant Module Trip 3.19
7147	7+5 ¹	Attendant Module Trip 3.19
7821	-	Thermal Feedback Active - 3.17
7825	-	Thermal Shutdown - 3.17
7902	-	Thermal Feedback Active - 3.17
All Other Codes	7 or 8	Possible Control System Trip - 3.7 & 3.8

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地址 - XXXXXXXXXXXX
XXXXXXXXXX
XXXXXXXXXX

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材積 - 40 X 25.5 X 15 (cm)

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控制器種類 (VSI, S-Drive, Egls, etc.)	PG 產品編號 (Dxxxxxx)	序號	退回原因
ECIS 110	D50238/1	25010184	10 BARS FLASHING
ECIS 110	D50238/1	45010195	WHILE CHARGING THE BATTERY, SCOOTER STILL CAN MOVE
TRUCHARGE	D50066	P5010968	LED CANNOT BE ILLUMINATED

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